



# Building monocyte derived dendritic cell behavior



## Experimental steps

This tab includes data related to wet-lab experiments, such as laboratory protocols and FASTQ files, done on monocyte derived dendritic cells to calibrate and validate intracellular models of monocyte derived dendritic cells.



## Bioinformatics steps

This tab includes data related to bioinformatics analysis, such as scripts and pipelines, used to analyze data gathered in the wet-lab part.



## Modeling steps

This tab includes data related to the intracellular mathematical model of dendritic cells, including literature used to build the models, links to scripts, results, etc.

# Experimental flow for coculture experiment

Monocyte derived dendritic cells

Sample metadata

PBMC Isolation Protocol

Differentiation and Protocol

GM-CSF+IL4

A. fumigatus

Transcriptome

Proteome

Collect at 0, 2,4,6,8 hours

RNA-SEQ

Protein measurement


Blot assays

ELISA

FASTQ files

Library Prep

Go to bioinformatics step



# Bioinformatics pipeline

Monocyte derived dendritic cells

Get Nextflow pipeline

FASTQ files for Dendritic Cells

MultiQC pretrim

Adapters FASTA file

Get reference genome

Analysis scripts

Pool samples from different runs

Cut adapters (Trimmomatic)

STAR aligner (index and align)

Qualimap, SAM tools  
MultiQC

DESEQ2, CHEA3,  
enrichment

Aligned reads

Results

Go to modeling step

# Mathematical modeling pipeline